



## Sonablate®500 world-wide unique

### FEATURES

and

### BENEFITS

Image-guided Therapy:  
Same confocal transducer for treatment  
and imaging

Precise planning and ongoing control  
during treatment; efficacy + safety  
Important safety and reference feature; you see  
what you do!

Multi-focal lengths probes

Treatment of ALL prostate zones within a  
40gr prostate without need to change the  
probe, precision and safety in dorsal parts!

High frequency of 4 MHz for treatment

Produces shorter and sharper lesions but with  
bigger diameter than lower MHz and less  
cavitation!

3D imaging planning

Allows better prostate visualization for prostate  
irregularities and ease of treatment planning.

"Simul-Plan" therapy planning (reference images)

Safety and feedback in treatment monitoring;  
allows doctor to view all treated tissue and  
compare with pretreatment tissue.

RIM (Reflectivity Index Measurement™)

Continuous Real-Time Rectal Wall Monitoring  
for safety – to assure protection of rectal  
wall and avoid damages/complications

Flexible power output/control  
Bypass mode

Power level changed/adjusted quickly according to  
actual needs/situations or completely suppressed

Flexible probe arm, adjustable  
to any surgical table

Precise and stable positioning of probe  
during treatment, high flexibility in adjusting

Lithotomy position for treatment

No special table/equipment necessary,  
whole equipment small with no big space  
requirement; flexibility in treating larger  
prostates and reaching anterior segments of  
prostate, no TUR necessary, patient comfort